

iOffice II™

A precision - engineered desk work progressive lens, **SWISSCOAT®iOfficeII™** is ideal for your desk work. It is an unique lens design which offers you more freedom of movement in the office. This design has extremely wide near, intermediate zone and clear vision field, they are the best choice for the wearer when working with computers, at work, at home and at play. This is a real individual near - vision lenses that guarantee you optimum comfort in wear.

- Ordering with progressive prescription
- Degression power automatically be calculated by software
- Choice for intermediate vision distance:
65cm ~ 200cm (5cm per step) on eye level

- * **Extended near and intermediate vision area**
- * **Customised with intermediate working distance**
- * **Available with Clear and Transitions Lens**

- Fitting Cross: center point of engraving axis
- Minimum vertical frame size: 26mm
- Minimum fitting height: 16mm
- 10mm from fitting cross to the top of inner frame
- Variable Inset Available

Progressive engraving:

Temple

ADD

2m


XX

Nasal >




Nasal Engraving XX / Supply Reading Power Range

Index \ Material		iReader II™ (1m)		iOffice II™ (2m)		iRoom II™ (4m)	
		Clear	Transitions® (Brown / Grey)	Clear	Transitions® (Brown / Grey)	Clear	Transitions® (Brown / Grey)
1.5 (1.498)	Nasal Engraving	E	T	E	T	E	T
	Order Code	1WR	1WRTB / 1WRTG	1WO	1WOTB / 1WOTG	1WM	1WMTB / 1WMTG
	Reading Power Range	+7.50D ~ -4.50D	+5.50D ~ -4.50D	+7.50D ~ -4.50D	+5.50D ~ -4.50D	+7.50D ~ -4.50D	+5.50D ~ -4.50D
1.58 Polycarbonate	Nasal Engraving	L	TL	L	TL	L	TL
	Order Code	1WRL	1WRTL B / 1WRTL G	1WOL	1WOTL B / 1WOTL G	1WML	1WMTL B / 1WMTL G
	Reading Power Range	+8.00D ~ 5.50D	+7.00D ~ -5.50D	+8.00D ~ 5.50D	+7.00D ~ -5.50D	+8.00D ~ 5.50D	+7.00D ~ -5.50D
1.53 Trivex	Nasal Engraving	V	TV	V	TV	V	TV
	Order Code	1WRVX	1WRTXB / 1WRTXG	1WOVX	1WOTXB / 1WOTXG	1WMVX	1WMTXB / 1WMTXG
	Reading Power Range	+6.00D ~ -5.50D	+6.00D ~ -4.50D	+6.00D ~ -5.50D	+6.00D ~ -4.50D	+6.00D ~ -5.50D	+6.00D ~ -4.50D
1.61	Nasal Engraving	H	TH	H	TH	H	TH
	Order Code	1WRH	1WRTHB / 1WRTHG	1WOH	1WOTHB / 1WOTHG	1WMH	1WMTHB / 1WMTHG
	Reading Power Range	+6.50D ~ -9.00D	+9.00D ~ -9.50D	+6.50D ~ -9.00D	+9.00D ~ -9.50D	+6.50D ~ -9.00D	+9.00D ~ -9.50D
1.67	Nasal Engraving	P	TP	P	TP	P	TP
	Order Code	1WRP	1WRTPB / 1WRTPG	1WOP	1WOTPB / 1WOTPG	1WMP	1WMTPB / 1WMTPG
	Reading Power Range	+11.50D ~ -8.50D	+9.00D ~ -9.50D	+11.50D ~ -8.50D	+9.00D ~ -9.50D	+11.50D ~ -8.50D	+9.00D ~ -9.50D
1.74	Nasal Engraving	U	TU	U	TU	U	TU
	Order Code	1WRU	1WRTUB / 1WRTUG	1WOU	1WOTUB / 1WOTUG	1WMU	1WMTUB / 1WMTUG
	Reading Power Range	+8.00D ~ -10.50D	+9.50D ~ -10.50D	+8.00D ~ -10.50D	+9.50D ~ -10.50D	+8.00D ~ -10.50D	+9.50D ~ -10.50D



Temple

LEFT LENS
(Concave)



Nasal

RIGHT LENS
(Concave)

Temple

