

iRoom II™

Whether at work, in transit or at home, typical modern city dwellers are spending most of their time indoors in front of their portable computing/telecommunications devices, an lifestyle evolution which creates special challenges and needs for eye care products.

SWISSCOAT® iRoom II™ is specifically designed to cater for this lifestyle evolution in the indoor habitat by offering an extremely comfortable indoor lens with tremendously extended middle vision zone to allow the wearer to work relaxedly in front of a computer while retaining perfect visual function in the typical indoor environment.

- Order data as general progressive lens by giving Far Vision power and Addition reading power
- The degressive power of iRoom II™ will be automatically calculated by iRoom II™ program
- The working zone of iRoom II™ – generally within 4 metres from the reading zone
(the actual distance of the working zone depends very much on the addition power and result of the iRoom II™ program calculation)
- “+” fitting cross located in between the two progressive engravings
- Variable Inset Available

iRoom II™ Progressive Marking:




iOffice II™



iReader II™




iRoom II™




Nasal Engraving XX / Supply Reading Power Range

Index \ Material		iReader II™ (1m)		iOffice II™ (2m)		iRoom II™ (4m)	
		Clear	Transitions® (Brown / Grey)	Clear	Transitions® (Brown / Grey)	Clear	Transitions® (Brown / Grey)
1.5 (1.498)	Nasal Engraving	E	T	E	T	E	T
	Order Code	1WR	1WRTB/ 1WRTG	1WO	1WOTB/ 1WOTG	1WM	1WMTB/ 1WMTG
	Reading Power Range	+7.50D ~ -4.50D	+5.50D ~ -4.50D	+7.50D ~ -4.50D	+5.50D ~ -4.50D	+7.50D ~ -4.50D	+5.50D ~ -4.50D
1.58 Polycarbonate	Nasal Engraving	L	TL	L	TL	L	TL
	Order Code	1WRL	1WRTLb / 1WRTLg	1WOL	1WOTLb / 1WOTLg	1WML	1WMTLb / 1WMTLg
	Reading Power Range	+8.00D ~ 5.50D	+7.00D ~ -5.50D	+8.00D ~ 5.50D	+7.00D ~ -5.50D	+8.00D ~ 5.50D	+7.00D ~ -5.50D
1.53 Trivex	Nasal Engraving	V	TV	V	TV	V	TV
	Order Code	1WRVX	1WRTXB / 1WRTXG	1WOVX	1WOTXB / 1WOTXG	1WMVX	1WMTXB / 1WMTXG
	Reading Power Range	+6.00D ~ -5.50D	+6.00D ~ -4.50D	+6.00D ~ -5.50D	+6.00D ~ -4.50D	+6.00D ~ -5.50D	+6.00D ~ -4.50D
1.61	Nasal Engraving	H	TH	H	TH	H	TH
	Order Code	1WRH	1WRTHB / 1WRTHG	1WOH	1WOTHB / 1WOTHG	1WMH	1WMTHB / 1WMTHG
	Reading Power Range	+6.50D ~ -9.00D	+9.00D ~ -9.50D	+6.50D ~ -9.00D	+9.00D ~ -9.50D	+6.50D ~ -9.00D	+9.00D ~ -9.50D
1.67	Nasal Engraving	P	TP	P	TP	P	TP
	Order Code	1WRP	1WRTPB / 1WRTPG	1WOP	1WRTPB / 1WRTPG	1WMP	1WMTPB / 1WMTPG
	Reading Power Range	+11.50D ~ -8.50D	+9.00D ~ -9.50D	+11.50D ~ -8.50D	+9.00D ~ -9.50D	+11.50D ~ -8.50D	+9.00D ~ -9.50D
1.74	Nasal Engraving	U	TU	U	TU	U	TU
	Order Code	1WRU	1WRTUB / 1WRTUG	1WOU	1WOTUB / 1WOTUG	1WMU	1WMTUB / 1WMTUG
	Reading Power Range	+8.00D ~ -10.50D	+9.50D ~ -10.50D	+8.00D ~ -10.50D	+9.50D ~ -10.50D	+8.00D ~ -10.50D	+9.50D ~ -10.50D



Temple

LEFT LENS
(Concave)



Nasal

RIGHT LENS
(Concave)

Temple

